

# Triamin Disinfection

Disinfection wipes for clinical use –  
alcohol free - no after rinse needed



Item no.: 83133

## Intended use

Triamin disinfection is intended to be used for disinfection of non-critical medical devices, which can be cleaned and disinfected at their point-of-use. Note that non-critical devices are those whose surfaces contact only intact skin or have no direct patient contact, yet may become contaminated during patient care.

## Contraindication

Do not use Triamin Disinfection for disinfection of semi-critical medical devices that contact mucous membrane or breached skin.

## Inherent safety

The pH-value and biocide concentration is lowered as much as possible in Triamin Disinfection to protect patient and users from extreme pH-values and hazardous residues. In addition, benzalkonium chloride and didecyldimethyl ammonium chloride are excluded from Triamin Disinfection because these biocides are sticky and difficult to remove from surfaces of medical devices.

## Protective measures during reprocessing

Triamin Disinfection ensures the disinfection of improperly cleaned devices by combining a high-performing viscose-cleaning substrate with a soil-tolerant biocide.

## Risk assessment

The risk assessment for Triamin Disinfection has been reviewed by Notified body Bureau Veritas Italy for compliance to the requirements of ISO EN 14971.

## Device classification

Class B  
Rule 15

## Material compatibility

Through the study of disinfectants in the literature and related in vitro data, we have found Triamin Disinfection to be compatible with hard surfaces made of metal, plastic, ceramics/glass, hard rubber and various proprietary plastic materials.



GDPMD  
CAB No: MDA/CAB 009

Certificate of Conformity: MYP120573  
MDA Reg.: GB9924620-39592

## Wet Wipe A/S

info@wetwipe.eu  
Tel: +45 70 266 244  
www.wetwipe.eu

CVR/VAT no: DK 30207726

## Aras Maju Hygiene

Level 7, Oasis Wing,  
Brunsfeld Oasis Tower 3  
No.2 Jalan PJU 1A/7A, Oasis Square,  
Oasis Damansara Petaling Jaya,  
47301 Selangor, Malaysia

Tel: +603-78485953  
E-mail: info@arasmaju.com

81-057, EN, Ver. 1  
Feb 2020

**wetwipe®**  
Developed in Denmark  
Made in Germany

# Documentation for Triamin Disinfection

## Device compatibility

Device specific compatibility documentation is generated through feedback from device manufacturers. Triamin Disinfection has already been entered in the positive list of GE Health Care. See: <https://www.gehealthcare.com/en/products/ultrasound/ultrasound-transducers>

## Biocompatibility

To demonstrate biocompatibility of our Triamin Disinfection wipes, a biological evaluation was performed by European Registered Toxicologist Dr. Anja Rämisch. Duration of skin contact limited to < 24h in the evaluation. No toxicological concern was identified. Therefore, no adverse effects are expected for the patient in contact to Triamin Disinfection residues via disinfected surfaces of medical devices.

## Clinical safety

To demonstrate clinical safety a clinical evaluation was performed by Dr. Iva Nestic and Dr. med. Isabelle Lang-Zwosta, according to MEDDEV 2.7.1 rev. 4.

## Isolation precautions

Some clinical situations may require the use of high-level disinfection (i.e., patients with norovirus infection, patients with Clostridium difficile infection).

## Shelf-life

The shelf-life for Triamin Disinfection is verified by testing the liquid extracted from real-time aged wipes 30 days after opening of packs.

## Ingredients

3-Aminopropyl dodecylamin, surfactants and foam-damper. Fragrance-free.

Color	Colourless
Odeur	Characteristic amine
pH-value	8-9

## Cloth material

High-performing viscose-cleaning substrate.

Biocidal efficacy	Exposure time
Bactericidal EN 13727, EN 16615	5 min
Yeasticidal EN13624, EN 16615	5 min
Vaccinia and blood-borne viruses DVV/RKI guideline	15 / 5 min

## Ordering information

Triamin Disinfection wipes are available for clinical use in one size:

	Size	Wipes/pack	Packs/carton	UMDNS-Code	PZN	Art.-Nr.
Mini	30 × 20 cm	25	20	15-281	13654602	83133



	Size	Pcs./package	Art.-Nr.
Wall mount mini	14,5 × 11,5 cm	1	30065

**Wet Wipe A/S**  
info@wetwipe.eu  
Tel: +45 70 266 244  
www.wetwipe.eu

CVR/VAT no: DK 30207726

**Aras Maju Hygiene**  
Level 7, Oasis Wing,  
Brunsfield Oasis Tower 3  
No.2 Jalan PJU 1A/7A, Oasis Square,  
Oasis Damansara Petaling Jaya,  
47301 Selangor, Malaysia

Tel: +603-78485953  
E-mail: info@arasmaju.com

**wetwipe®**  
Developed in Denmark  
Made in Germany